

Whatever.....your shape.....size or type.....of stoma....we can help....you to find.....the perfect fit.

HARMONY DUO with Flexifit®			
STANDARD FLANGES with Flexifit and Aloe Vera extracts			
SIZE	ORDER CODES	PACK SIZE	
13-32mm	FHD1332	10	
13-50mm	FHD1350	10	
13-70mm	FHD1370	10	
FLEXIFIT FLANGES with hydrocolloid surround and Aloe Vera extracts			
SIZE	ORDER CODES	PACK SIZE	
13-32mm	FHDF1332	10	
13-50mm	FHDF1350	10	
13-70mm	FHDF1370	10	
LARGE CLOSED POUCHES (Colostomy)			
SIZE	ORDER CODES Opaque	PACK SIZE	
13-32mm	HDCL1332	30	
13-50mm	HDCL1350	30	
13-70mm	HDCL1370	30	
STANDARD CLOSED POUCHES (Colostomy)			
SIZE	ORDER CODES Opaque	PACK SIZE	
13-32mm	HDC1332	30	
13-50mm	HDC1350	30	
13-70mm	HDC1370	30	
MINI CLOSED POUCHES (Colostomy)			
SIZE	ORDER CODES Opaque	PACK SIZE	
13-32mm	HDCM1332	30	
LARGE DRAINABLE POUCHES (Ileostomy)			
SIZE	Transparent	ORDER CODES Trans/Overlap	PACK SIZE
13-32mm	-	HDDL1332	30
13-50mm	HDDL1350*	HDDL1350	30
13-70mm	HDDL1370*	HDDL1370	30
STANDARD DRAINABLE POUCHES (Ileostomy)			
SIZE	ORDER CODES Trans/Overlap	PACK SIZE	
13-32mm	HDD1332	30	
13-50mm	HDD1350	30	
13-70mm	HDD1370	30	
SMALL DRAINABLE POUCHES (Ileostomy)			
SIZE	ORDER CODES Trans/Overlap	PACK SIZE	
13-32mm	HDDS1332	30	
13-50mm	HDDS1350	30	
13-70mm	HDDS1370	30	

\*Post-op pouches with NO filter



# We've found a perfect fit for every body.

**Confidence<sup>®</sup> Natural**  
WITH ALOE EXTRACTS

If you have a hernia or uneven skin around your stoma, it can be difficult to find a pouch that fits. Available in a choice of drainable, closed or urostomy pouches, Confidence<sup>®</sup> Natural features our unique five-sided Flexifit<sup>®</sup> wafer, which is specially shaped to fit more closely with the body's natural contours, helping to reduce leaks and irritation. The hydrocolloid is even infused with extracts of Aloe to help keep skin smooth and soft. The perfect fit never felt so good.

### Natural softness


Skin-friendly hydrocolloid infused with Aloe Vera for softer skin

### Natural protection

Skin-friendly hydrocolloid, research into which is validated by the British Skin Foundation

### Natural security (urostomy only)

Special urostomy adhesive keeps pouch firmly in place and breaks down more slowly, which could provide a longer wear time



DRAINABLE POUCH FOR ILEOSTOMISTS



FLEXIFIT WAFER ON A HERNIA



CLOSED POUCH FOR COLOSTOMISTS



UROSTOMY POUCH WITH TAP

CONFIDENCE NATURAL with Aloe extracts				
LARGE CLOSED POUCHES (Colostomy)				
SIZE	Transparent	ORDER CODES	Trans/Overlap	PACK SIZE
13mm	NLT13*		NL13*	30
STANDARD CLOSED POUCHES (Colostomy)				
SIZE	Transparent	ORDER CODES	Trans/Overlap	PACK SIZE
13mm	NT13*		N13*	30
25mm	-		N25	30
28mm	-		N28	30
32mm	-		N32	30
35mm	-		N35	30
38mm	-		N38	30
41mm	-		N41	30
MINI CLOSED POUCHES (Colostomy)				
SIZE	Transparent	ORDER CODES	Trans/Overlap	PACK SIZE
13mm	-		NM13*	30
LARGE DRAINABLE POUCHES (Ileostomy)				
SIZE	Transparent	ORDER CODES	Trans/Overlap	PACK SIZE
13mm	NDL13*		NDL13*	30
STANDARD DRAINABLE POUCHES (Ileostomy)				
SIZE	Transparent	ORDER CODES	Trans/Overlap	PACK SIZE
13mm	NDT13*		ND13*	30
25mm	-		ND25	30
28mm	-		ND28	30
32mm	-		ND32	30
35mm	-		ND35	30
SMALL DRAINABLE POUCHES (Ileostomy)				
SIZE	Transparent	ORDER CODES	Trans/Overlap	PACK SIZE
13mm	-		NDS13*	30
STANDARD POUCHES (Urostomy)				
SIZE	ORDER CODES	Trans/Overlap	PACK SIZE	
13mm*	NU13			10
25mm	NU25			10
28mm	NU28			10
32mm	NU32			10
35mm	NU35			10
SMALL POUCHES (Urostomy)				
SIZE	ORDER CODES	Trans/Overlap	PACK SIZE	
13mm	NUS13*			10

\*Starter hole wafers can be cut up to 70mm wide and 55mm deep


**Harmony<sup>®</sup> Duo**  
with Flexifit<sup>®</sup> Flange & Aloe Extracts



STANDARD FLEXIFIT FLANGE



FLEXIFIT FLANGE WITH HYDROCOLLOID SURROUND



LOCATOR



FLEXIFIT FLANGE ON A HERNIA

Those who enjoy the convenience of a two-piece system can also benefit from an improved fit, thanks to Harmony<sup>®</sup> Duo. These discreet pouches now feature the unique Flexifit flange, which fits closely to any body shape; and a softer hydrocolloid surround, giving you greater freedom of movement. With stronger adhesive to keep your pouch firmly in place and Aloe extracts to help soothe the skin, Harmony Duo offers our most comfortable and secure fit ever.



DRAINABLE POUCH FOR ILEOSTOMISTS

### NEW Flexifit wafer

Unique, five-sided Flexifit wafer fits closely to any body shape.

### NEW Stronger adhesive

Stronger and stickier adhesive for extra security and peace of mind.

### NEW Aloe extracts

Skin-friendly hydrocolloid infused with extracts of Aloe Vera for softer skin.

### Easy application

A locator on the wafer matches up with a notch on the pouch for simple, accurate positioning every time.

### Shaped to fit

Excellent flexibility and adhesion, ideal for hernias or uneven skin.



CLOSED POUCH FOR COLOSTOMISTS