

SALTS ADHESIVE REMOVER

DESCRIPTION	ORDER CODES	SIZE
Spray	WAP	50ml
Wipes	WA1	30 sachets



Remove and refresh

Dual action adhesive remover

For free samples, please contact:



TOLL FREE 1300 784 737 (Australia)



@SaltsHealthcare

TOLL FREE 0800 100 146 (New Zealand)



SaltsHealthcare



ostomy@ainscorp.com.au



www.ainscorp.com.au

CE ©Registered trade marks of Salts Healthcare Ltd. ©Salts Healthcare Ltd 2018. Products and trade marks of Salts Healthcare Ltd are protected by UK and foreign patents, registered designs and trade marks. For further details, please visit www.salts.co.uk

RM764960 03/18

DETAILS REQUIRED

salts
HEALTHCARE



TWO KEY BENEFITS IN ONE PRODUCT— SALTS ADHESIVE REMOVER REMOVES AND REFRESHES

As well as gently and effectively removing residue, Salts Adhesive Remover Spray has a no-sting formula with a natural, peppermint fragrance which helps to mask odour. So you can feel fresh after every bag change.

Salts Adhesive Remover Spray also features 'bag-on-valve' aerosol technology, which means it can be sprayed in any direction, even upside down, helping you to reach the most awkward places.

- ▶ No-sting formula
- ▶ Bag-on-valve aerosol technology means it sprays at all angles – 360°
- ▶ Reduced wastage – up to 100% product emptying
- ▶ Quiet spray for discreet use
- ▶ Won't feel cold on skin
- ▶ Fresh peppermint fragrance, helps mask odour – 94% of patients agreed*



SALTS ADHESIVE REMOVER WIPES

Now our best-selling no-sting adhesive remover formulation with natural peppermint fragrance is available in single-use wipes. Our new wipes are available in boxes of 30, and a convenient tear notch on each sachet makes them easy to use.



* Product assessment. Data on file at Salts Healthcare Ltd.

THE BENEFITS OF BAG-ON-VALVE TECHNOLOGY

THE ACTUATOR

Sprays at any angle (360°), no pumping required

THE VALVE

Ensures controlled flow and consistent spray

THE CANISTER

Compact and convenient 50ml aluminium can

THE BAG

Separates product from propellant for longer shelf life

THE PROPELLANT

Squeezes product out using compressed air/nitrogen

THE PRODUCT

Less-chilling product discharge – ideal for sensitive skin



98%

said the adhesive remover **could be sprayed upside down***

96%

said **bag application was easy** after using the spray*

84%

said that the spray was the same as or **better than their usual product***