

SALTS FLANGE EXTENDER WITH ALOE



PICTURE FRAMING FOR EXTRA SECURITY

Features:

- ▶ Helps secure the edges of the wafer, providing added security when needed.
- ▶ Aloe Vera extracts to help soothe and protect
- ▶ One size fits all
- ▶ 120° – use as many strips as required (3 strips will form a full circle)
- ▶ Extra security with any type of bag
- ▶ Non-allergenic – great for people who can't use microporous tapes



HOW TO APPLY



1 Remove the film backing.



2 Position as required and smooth onto the skin.

SALTS FLANGE EXTENDER WITH ALOE

SIZE	ORDER CODE	PACK SIZE
One size fits all	SPHA2	30